

Vibrating Needle During Venipuncture Reduces Insertion Force and Yields Lower and Less Variable Average Corticosterone Levels in Rodents

R.S. Clement¹, E.L. Unger², S.A. Cavigelli³, R.M. Sheehan¹,
R.B. Bagwell¹, V.A. Kellogg¹, M.L. Mulvihill¹

¹ Actuated Medical, Inc.
310 Rolling Ridge Dr.
Bellefonte, PA 16823

² Department of Biology
Lebanon Valley College
Annville, PA 17003

³ Department of Biobehavioral Health
The Pennsylvania State University
University Park, PA 16802

Pat. actuatedmedical.com/ip



Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Support

This research is supported by Small Business Innovation Research (SBIR) Phase I and II Grants from the National Institutes of Health, National Institute on Aging.

“The project described was supported by Award Number R44AG037214 from the National Institute on Aging. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institute on Aging or the National Institutes of Health.”

Pat. actuatedmedical.com/ip



Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Animal Research Needs

+ Animal Research Needs

- + Described in the 2009 SBIR Omnibus Solicitation of the National Institutes of Health, Division of Aging Biology.
 - + *Development of minimally-perturbing techniques for collecting blood from mice, rats, and other animals several times a day in sufficient quantities for measurement of hormone levels and other circulating factors in you*
- + The 3 R's: Reduction, Refinement, Replacement
 - + *Less stressful sampling (Refinement)*
 - + *Less variability = less animals (Reduction)*

Pat. actuatedmedical.com/ip




Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Basis of Approach

- + Premise of the Solution: Vibration.
 - + Gate Theory of Pain – describes anesthetic effect of vibration.
 - + Mosquitos – drawing blood in nature.



Pat. actuatedmedical.com/ip

Proposed Solution

Vibrating Needle for Venipuncture, GentleSharp



Phase I SBIR Hypothesis:

- + An actuated resonance-assisted lancet introduction device will significantly reduce insertion force (>50%), leading to less stressful blood sampling in rodents, without causing additional tissue damage.
- + The hypothesis was tested in a serial blood sampling study with Sprague Dawley Rats.

Pat. actuatedmedical.com/ip



Innovative motion + Positive outcomes


©2013 Actuated Medical, Inc.

In vitro Testing, Cadaver Rat Tails

Insertion Force Measurement Set-up

Displacement
Sensor Conditioner

Data Acquisition
Hardware




Current
Probe

Control
Unit

Force Gauge

Sample
Tray

Prototype



Needle Attached
to Prototype

Force
Gauge

Rat Tail
Segment

Pat. actuatedmedical.com/ip




Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

In vitro Insertion Force vs. Time Plots

Vibrated Needle Exhibits Lower Insertion Force



Pat. actuatedmedical.com/ip





Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

In vitro Peak Insertion Force Comparison

Vibrated Needle Exhibits Lower Peak Force



+ Testing performed on 6 rat tails at three locations.

- + N = 18 on/18 off.
- + Error bars = standard deviation.

Pat. actuatedmedical.com/ip



Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

In vivo Serial Blood Sampling Experiment

Sprague Dawley Rat Study, Tail Site



- + Research was conducted under an IACUC approved protocol in the Department of Nutritional Sciences, the Pennsylvania State University.
- + PI: Dr. E. Unger (Phase I SBIR).

- + Groups
 - + Treatment (on): 10 subjects.
 - + Control (off): 9 subjects.
- + Protocol
 - + Sample days occurred 3x at 1 week intervals.
 - + 3 blood samples attempted in each rat on each sample day (each sample separated by ~1 hr).
- + Data
 - + Corticosterone concentration.
 - + Number of attempts required for each sample (success/failure).
 - + Vocalization/Movement (Likert).
 - + Presence/absence of hematoma.

Pat. actuatedmedical.com/ip




Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Stress Hormone Concentration Comparison

Subjects with Vibrated Needle Exhibit Lower Cort. Levels



Difference: **9.7%** **49.2%** **65.2%**

Pat. actuatedmedical.com/ip




Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Over Study Duration Stress Hormone Comparison

Subjects with Vibrated Needle Exhibit Less Cort. Variation

Average Shift from Initial Baseline (by Subject) – Cort.



Pat. actuatedmedical.com/ip




Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Plasma Stress Hormone Comparison

Subjects with Vibrated Needle Exhibit Lower Cort. Variability



Pat. actuatedmedical.com/ip




Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Cage Mate Comparison

Subjects with Vibrated Needle Exhibit 33.7% Lower Cort.



- + Pair-housed rat comparison.
 - + One sampled always with device on and the other with device off.
 - + Cage 1: both subjects had device on; Cage 3: both subjects had device off, thus not shown.
- + In every case, the rat sampled with device on had lower average corticosterone release.
- + By a cage by cage basis, the rat sampled with device on had **33.7% lower corticosterone levels** as compared to off.

Pat. actuatedmedical.com/ip




Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Behavioral Comparison

Subjects with Vibrated Needle Exhibit Lower Vocalization & Movement



Pat. actuatedmedical.com/ip




Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Focus Group Survey Results

End-users Agreed Reducing Subject Stress Important



- + Veterinarians, researchers, and technicians (N=9) participated in focus groups.
- + The participants inserted needles into several models, including a cadaver rat tail, with and without Vibrating Needle on.
- + The participants then answered a questionnaire to evaluate their experience with the Vibrating Needle, GentleSharp, and its potential.

Pat. actuatedmedical.com/ip



Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Study Summary

- + Vibration reduced insertion force into rat tails by up to 68%.
- + Up to **65% less stress hormone** ($p < 0.01$) observed in the blood samples obtained with the prototype turned on.
- + Up to **40% less average standard deviation of individual corticosterone** measured in samples with vibration compared to without vibration.
- + Statistically significant reduction in behavioral markers of stress (vocalization and movement) with vibration during venipuncture.
- + Focus Group of potential end-users all agreed very strongly that the vibrated needles were easier to insert.
- + Focus Group of potential end-users agreed that reducing animal discomfort and stress during blood sampling was important and that after reviewing the data and trying device believed it has value.
- + Supports Reduction and Refinement of animal studies that require blood sampling.

Pat. actuatedmedical.com/ip



Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Current/Future Work

- + Ongoing study at Penn State University with mice to evaluate other sampling sites.
- + Ongoing study at Medical University of South Carolina to evaluate tail blood sampling of mice.
- + 12-month longitudinal serial blood sampling study at the Pennsylvania State University, starting December 2013.
- + The Vibrating Needle Device, GentleSharp anticipated product launch January 2014.

Pat. actuatedmedical.com/ip



Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.

Thank You to:

- + Nancy Nadon, Ph.D.,
Program Officer, Biological
Resources Program, National
Institute on Aging.
- + National Institutes of
Health, National Institute
on Aging SBIR Program.
- + And to you for your time
today!

Questions??

 GentleSharp™



GentleSharp can be
previewed at **Booth 2143**
in the Exhibit Hall.

Pat. actuatedmedical.com/ip



Innovative motion + Positive outcomes

©2013 Actuated Medical, Inc.